# neno<sup>®</sup>

## bene

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PL

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**Producent:**  
Manufacturer's data:  
Vapo Healthcare Co., Ltd.  
Southern unit of third floor, building B, No.99  
Yuda West Rd, High tech district, KunShan,  
Suzhou, Jiangsu Province, China.

**Autoryzowany przedstawiciel:**  
EU authorized representative:  
Share Info GmbH  
Heerdter Lohweg 83, 40549 Düsseldorf

# nenobene

## bene

USER'S MANUAL

Dear Customer,  
thank you for buying our product, the Neno Bene mesh nebulizer.  
Please read the user manual before use.

### 01. DESCRIPTION

**Method of operation:** Thanks to the self-exciting oscillation circuit, the ceramic walls together with a metal mesh and a built-in nebulizing head guarantee successful oscillation under a high frequency. The device structure ensures smooth flow of the medicine solution into the micromesh of the metal mesh and nebulization into an aerosol form, which flows into patient's respiratory system through the inhalation mask to successfully inhale medicine.

#### Symbols:

Please read the information on the signs listed in the table to avoid user mistakes and to become familiar with the safety measures.

The device should be utilized in designated municipal waste areas

Class 2


User manual

BF type parts

**SN** Serial number given by the manufacturer

#### Accessories

Before installation and configuration of the device, check if the device is in good shape and if the quality and quantity of the accessories matches the quantity listed in the user manual. Check, if there are any device defects present, if so, do not use the device, but contact the manufacturer or seller instead.

1. Medicine Cup
2. Mask
3. String
4. Main Unit
5. Main Unit Power Button
6. Main Unit Connection Cable
7. Type-C Adapter Head

### 02. SPECIFICATION

|  |  |
|--|--|
| <b>Power :</b>                           | DC 5V 1A   |
| <b>Aerosol efficiency:</b>               | ≥0.20 ml/min   |
| <b>Particle size:</b>                    | MMAD approx. 5µm   |
| <b>Sound level:</b>                      | ≤35dB  |
| <b>Ultrasonic frequency:</b>             | Approx. 114KHz   |
| <b>Medicine cup capacity:</b>            | 8ml  |
| <b>Medicine cup temperature range:</b>   | ≤60°C  |
| <b>Medicine cup weight:</b>              | Approx. 11g  |
| <b>Dimensions:</b>                       | 42mm x 34mm x 40mm   |
| <b>Safety measures :</b>                 | Do not use the device near flammable, intoxicating gases mixed with air, or flammable, intoxicating gasses mixed with oxygen or nitric oxide |
| <b>Application method:</b>               | Medicine in nebulized form   |
| <b>Contraindications:</b>                | None   |
| <b>Electrical safety requirements:</b>   | Internal power mechanisms, type BF application parts   |
| <b>Work mode:</b>                        | Continuous   |
| <b>Work environment:</b>                 | Temperature: +10~+40°C<br>Humidity: 30~85%RH<br>Pressure: 86~106 kPa   |
| <b>Storage and transport conditions:</b> | Temperature: -20~+50 C<br>Humidity: 30~85%RH<br>Pressure: 86~106 kPa   |

### 03. SAFETY

**WARNING:**  
-When using a certain kind, dose, and pattern of medicine intake, make sure to adhere to doctor's recommendations. If failed to adhere to the doctor's recommendations, your symptoms may worsen.  
-Clean and disinfect the medicine cup and inhalation mask after every use.  
-Clean and disinfect device parts and then carefully dry them, and store in a dry, clean place to prevent bacterial growth.  
-Nebulizer is designed for individual use. It is not recommended to use the device by more than one person.  
-Keep the device out of reach of children and infants to avoid choking hazard. If swallowed, consult your doctor immediately.  
-Do not wash or immerse the device under water, especially hot. Risk of electrical installation damage, short circuit, or electrocution. Direct water contact can make the device useless.  
-Neno Bene is a professional nebulizer. Adults and children who are unable to use the device on their own should use the device under a guardian supervision.


#### ATTENTION:

-Do not inhale clear water from the medicine cup. Your symptoms may worsen.  
-Do not drop medicine on the main unit. If you drop the medicine, immediately wipe it with a gauze.  
-Using the main unit when wet can cause problems or even electrical shock.  
-Do not pierce or pick the metal mesh with a q-tip or a hair pin, it can cause damage to the set and make it unusable.  
-Do not drop or press forcefully on the medicine cup. It can cause damage or even electrical shock.  
-To ensure proper use of the device, read the manual carefully and adhere to the instructions. If you have any questions, consult the supplier or the manufacturer.  
-If the storage and transport temperature is lower than 5°C, the device should return to room temperature before use.  
-Please read the package leaflet before use to ensure that the medicinal product is suitable for the mesh nebulizer.

### 04. USE GUIDE

#### A. CONNECT THE POWER CABLE:

- The device uses an external power system (5V 1A). The main unit has an USB input and a power system (5V 1A). Insert the main unit USB plug into the power input. If you have a type C adapter head purchased, you can use it with the main unit. The main unit uses an external power supply with battery case.  
**Attention! Rechargeable batteries cannot be used in the battery case!**


#### B. BATTERY LIFE AND REPLACEMENT:

**1. When using alkaline batteries:**  
- The unit can be used for approximately 3 days (if used 30 minutes a day)  
- When the indicator light changes to red, the batteries are running low. Replace the batteries.  
- If you are planning not to use the device for a long time (3 months or more) take the batteries out.  
- If you leave the batteries inside the device, there is a risk of battery leakage and device damage.  
- Do not insert the batteries inside the device using the wrong polarities.


#### C. FILL THE MEDICINE CUP:

Open the medicine cup, fill it with medicine and close the lid tightly. Before using high density or strong medicine consult your doctor.

#### D. CONNECT THE MEDICINE CUP:

Insert the other end of the power cable into the port in the medicine cup, and make sure that it is mounted properly. Do not abruptly shake the medicine cup when filled with liquid.

#### E. MOUNT THE INHALATION MASK:


#### F. PUT THE MASK ON:

#### G. START THE INHALATION:

- If the battery case is used to power the device, make sure that it is turned on and that the indicator light is on. If an internal power supply is used, make sure, that the power supply is turned on and that it is compatible with the output voltage 5V 1A. Press the power button to start the inhalation.


-Take a deep breath to properly inhale the medicine.

-When using high density medicine, the nebulization can take longer than usual.  
-The sound coming out of metal sheet during the nebulization process is a normal phenomenon.  
-If the liquid has high surface activity, a small quantity of the mixture can flow out of the metal sheet.  
-If the medicine cup runs out of liquid, press the button on the device to avoid potential damage of the metal mesh.  
-If the medicine cup is running out of medicine, it is recommended to slightly tilt the cup, so the liquid will flow towards the metal mesh.  
-During mask inhalation, do not cover the air vents of the mask with your hands or other objects.

-During aerosol inhalation, do not cover the air vents on the lid with your hands or other objects to not influence the aerosol stream.  
-If the inhalation is causing you discomfort, immediately turn off the device and consult with your doctor.

#### H. STOP THE INHALATION


Press the power button on the device to stop inhalation. If the power source is the battery case, press the power button to turn off the battery case. Wait until the indicator lamp will turn off and disconnect the main unit from the battery case. If the device is powered by an external power source, turn it off, and unplug the main unit from the external power source.

### 05. CLEANING AND DISINFECTION

#### CLEANING THE MASK AND MEDICINE CUP

1. Take out the inhalation mask and the medicine cup out of the device.
2. Open the medicine cup and pour out the contents.
3. Add a small amount (2-5ml) of water into the medicine cup and cover it with the lid. Shake the cup lightly so that the medicine dissolves fully in the water. Open the medicine cup and pour out the contents.
4. Wash the mask with water.
5. After washing, dry the mask thoroughly.
6. Store all the nebulizer parts in a dry, clean place.

#### DISINFECTING THE MASK AND MEDICINE CUP:

1. Immerse the medicine cup in medical alcohol (75%) for no less than 3 minutes, then take the cup out of the alcohol and dry it thoroughly. Store the cup in a clean place.
2. Immerse the mask in medical alcohol (75%) for no less than 3 minutes. After cleaning, dry the mask thoroughly and store it in a clean place.
3. After disinfecting, rinse the device with water to avoid accumulation of medical alcohol residue. Make sure to disinfect the medicine cup after each use. In case of permanent staining change the device parts.

#### MASK AND MEDICINE CUP REPLACEMENT:

Medicine cup should be replaced after around 20 hours of use. In any other case your health and safety, functional capacity and safe inhalation will not be guaranteed. If the medicine cup is damaged or is working improperly, it must be replaced.

### 06. MAINTENANCE

- Check the device before every use. Do not use the nebulizer if it is damaged. Using the device when damaged can cause harm or danger.
- Store the device in a cold, dry and ventilated room. Keep away from fire and hot surfaces to avoid risk of battery explosion.
- If you have any problems with the device, like function settings, maintenance or using the device, contact the manufacturer. Do not deconstruct or fix the device on your own.
- In case of any occurring problems, contact the manufacturer.

### 07. TROUBLESHOOTING

| Problem  | Possible cause                              | How to fix  |
|--|---|---|
| Nebulization power is too weak                           | Membrane is dirty                           | Clean the medicine container. If the problem can not be fixed by the method above, replace the container with a new one.            |
| The device is not nebulizing                             | The main unit is not connected properly     | Unplug the main unit from the power source and connect it back again.   |
| The main unit light is on, but the device is not working | The main unit cable is not working properly | Replace the main unit with a new one.   |
|  | Device membrane is very dirty               | Clean the medicine cup according to the user manual. If there are problems after cleaning, replace the medicine cup with a new one. |
| Fluid leak   | Medicine cup is damaged                     | Replace the medicine cup with a new one.  |

### 08. DISPOSAL

DO NOT throw out the batteries with municipal waste. Throw out the batteries according to the local law on waste utilization (for example designated waste collection points). The device is made from plastic and stainless steel. Make sure to utilize materials according to the binding waste utilization law.

### 09. ELECTROMAGNETIC COMPATIBILITY

#### Guidance and manufacturer's declaration

1. Device parts should be installed and aligned according to the information included in the ATTACHED DOCUMENTS.  
2. IMMUNITY LEVEL TESTS for elementary rules of working of medical equipment and their software should be picked based on a big probability of maintaining security and overall productivity, and according to recommendations of professional health facilities, house environment and other special designated places of use.  
3. HOUSE ENVIRONMENT is a place in which the patient is living, or a place in which the patient is residing, except professional care facilities with trained medical professionals.  
FOR EXAMPLE: As indicated in table 6 IEC 60601-1-2: 2014 for medical equipment, a standard cellphone with a maximum output voltage of 2W stops at d=3.3m with an immunity level of 3V/m.

#### A1. Electromagnetic emissions for all devices and systems

#### Guidance and manufacturer's declaration - Electromagnetic emissions

The device is designed to use in electromagnetic environments listed below; patient using the nebulizer should make sure that they are in the correct environment.

| Emission test                    | Realized standard | Electromagnetic environment - tips  |
|----------------------------------|-------------------|---|
| RF CISPR 11 Emissions            | Group 1           | Nebulizer uses RF energy only for internal functions. Its RF emissions are very low and do not disrupt other devices from working.                            |
| RF CISPR 11 Emissions            | Class B           | The device is optimal to use at any place different than home and public environments, with low power or a mains power supply built in a domestic house grid. |
| Harmonic emissions IEC 61000-3-2 | Class A           |   |
| Harmonic emissions IEC 61000-3-2 | Compatibility     |   |

#### A2. Electromagnetic immunity - household environments.

#### ACCESSORIES AND SYSTEMS

#### Guidance and manufacturer's declaration - Electromagnetic immunity

The device is designed to use in an electromagnetic environment listed below; patient using the nebulizer should make sure that they are in the correct environment.

| Immunity test  | Measured level IEC 60601                            | Compatibility level                                 | Electromagnetic environment - tips  |
|--|---|---|---|
| Electrostatic emissions (ESD) IEC 61000-4-2            | ±8 kV contact<br>±2 kV, ±4 kV,<br>±8 kV, ±15 kV air | ±8 kV contact<br>±2 kV, ±4 kV,<br>±8 kV, ±15 kV air | The floors should be wooden, concrete or ceramic. If the floor is covered with a synthetic material, the relative humidity should be at least 30% |
| Electrical fast transient/burst 61000-4-4              | ±2kV for power lines                                | ±2kV for power lines                                | Power quality network It should be like this same as typical environment hospital or commercial.  |
| Surge IEC 61000-4-5                                    | ±1kV line to line<br>±2kV line to ground            | ±1kV line to line<br>±2kV line to ground            | Mains power quality should be the same as in a standard hospital or commercial environment.   |
| Power frequency (50/60Hz) magnetic field IEC 61000-4-8 | 30A/m   | 30A/m   | Power frequency of magnetic fields should be at levels characteristic for typical hospital and commercial environments.                           |

| Immunity test  | Measured level IEC 60601  | Compatibility level   | Electromagnetic environment - tips   |
|--|---|---|--|
| Voltage dips, short interruptions, and voltage variations on power supply input lines IEC 61000-4-11 | 0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° | 0% UT; 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° | Quality of the mains power supply should be typical for a hospital and commercial environments. If the user requires continued operation during power outages, it is recommended to power the device from an uninterrupter power supply.   |
| Conducted RF IEC 61000-4-6   | 3 Vrms 150Kz to 80MHz   | 3 Vrms  | The manufacturer should consider reducing the minimum separation distance based on risk management, and using higher immunity test levels, that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher immunity test levels shall be calculated using the following equation:<br>$E = \frac{B}{d} \sqrt{P}$ Where: "P" is the maximum power in watts, "d" is the minimum separation distance in meters and "E" is the immunity test level in V/m |
| Radiated RF IEC 61000-4-3  | 10 V/m 80 MHz to 2,7 GHz  | 10V/m   |  |

ATTENTION: Check the mains power supply before using the device

### 10. WARRANTY

The product is covered by a 24-month warranty. The terms of the guarantee can be found at: <https://nenop.pl/gwarancja>

Complaints should be reported using the complaint form located at: <https://nenop.pl/kontakt>

Specifications and contents of the kit are subject to change without notice. We apologise for any inconvenience.

KGK Trade declares that this device Neno Bene is compliant with significant requirements of directive 2014/53/EU. The text of this declaration may be found on the website: <https://nenop.pl/download/DOC/deklaracja-CE-Neno-Bene.pdf>

Last update of the user's manual: 01.02.2021

**Importer:**  
KGK Trade sp. z o.o. sp. k.  
Ujastek 5b  
31-752 Cracow  
Poland

Made in P.R.C.

The crossed out trash can symbol indicates that unusable electrical or electronic devices, its accessories (such as power supplies, cords) or components (for example batteries, if included) cannot be disposed of alongside with household waste. In order to dispose of the devices or its components (for example, batteries) deliver the device to the collection point, where it will be accepted free of charge. Disposal is subject to the recast version of the WEEE Directive (2012/19 / EU) and the Directive on batteries and accumulators (2006/66 / EC). Proper disposal of the device prevents degradation of the natural environment. Information about the collection points of the facilities is issued by the competent local authorities. Incorrect disposal of waste is subject to penalties provided for by the law in force in the given area.

**Manufacturer's data:**  
Vapo Healthcare Co., Ltd.  
Southern unit of third floor, building B, No.99  
Yudai West Rd. High tech district, KunShan,  
Suzhou, Jiangsu Province, China.

**EU Authorized Representative:**  
Share Info GmbH  
Heerdter Lohweg 83, 40549 Düsseldorf