# neno<sup>®</sup>

## medic T05

YK-IRT2

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**Wytwórca:**  
  
Xuhou Yongkang Electronic Science Technology Co., Ltd.  
4F building C8, 40, Jingshan Road, Economic and  
Technological Development Zone,  
221000 Xuzhou, PEOPLES REPUBLIC OF China

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BY **XBlitz**

**neno<sup>®</sup>**  
www.neno.pl  
Manufactured in China  
Designed by Harens

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YK-IRT2

USER'S MANUAL

**Dear Customer,**  
thank you for choosing Neno Medic T05. The device you purchased is an infrared thermometer that measures your body temperature or object using an infrared light wave sensor.  
**Before using the device, please read the instructions below.**

**5. View measurements stored in memory**  
When the device is turned on, press the "UP" button or the "DOWN" button to view the measurements stored in the memory. Press "DOWN" to view records from the oldest and "UP" to view records from the last measurement taken. The higher the number of the recording, the earlier the measurement. In other words, measurement 01 is always the last measurement taken by the device. If there is no stored measurement in the device memory, the measurement number will be displayed normally, but instead of showing measured temperature, the display will show ---. The thermometer can memorize up to 34 temperature measurements. If you take more measurements, the oldest saved result will be deleted to make room in memory for the new measurement.

**6. Clear measurement memory**  
When the device is turned off, hold down the MODE button and press the device trigger. "CLR" should appear on the display, indicating that all records in the device's memory have been deleted. The device will turn off two seconds after the message appears on the display screen.

**7. Change the temperature unit**  
Turn on the device. Hold down the MODE button until "F1" appears on the display. Press the "DOWN" button to switch the measurement unit between "C" and "F".

**Warning!** After pressing the MODE button twice, the message "F2" will show on the display screen. The "F2" message means that the service mode is turned on. The value displayed with the "F2" message should be 0.0. Changing this value when service mode is turned on will result in voiding the warranty!

### 5. DEVICE ERRORS

Message	Meaning	Possible reason and solution
Hi	Measured temperature too high	Body temperature measurement mode: The temperature of the test person is above 43 °C
Lo	Measured temperature too low	Body temperature measurement mode: The temperature of the test person is below 34 °C. Make sure you follow the device's instructions for use to ensure accurate measurement.
(blinks)	Low battery	Batteries in the device close to discharge. Replace the batteries with new ones.
	The thermometer turns off automatically after a moment of inactivity, this is the intended operation of the device	Restart the device by pressing the trigger
Blank screen	Batteries installed wrong	Make sure the battery poles are facing in the correct direction
	Discharged batteries	Replace batteries with new ones
	The display still does not start despite following the instructions above	Contact your device's service or seller

### 01. SAFETY MEASURES

- Do not use the infrared thermometer in ways other than the ones described in the manual. The thermometer is suitable both for home use and as a medical device.
- Do not immerse the thermometer in water or other liquids. When cleaning the device, refer to the instructions in the CLEANING AND STORAGE section.
- The thermometer should be stored in a dry, clean place, away from sunlight. The thermometer works best at a temperature range of 15-40°C and with air humidity in the range of 30-85%RH.
- Do not touch the thermometer sensor with your fingers.
- Obstacles such as sweat, hair hats or bandanas can lower the measured temperature. Make sure that nothing is blocking the path between the sensor and bare skin of the test subject.
- Do not drop the product, disassemble the device, or make repairs and modifications on your own.
- Keep the thermometer away from strong electrostatic fields or magnetic fields to avoid incorrect measurements.
- If you experience any problems, stop using the device and contact the seller or service.
- Do not dispose of the product or battery in the mixed municipal waste container. Comply with the laws regarding the disposal of electronic equipment and batteries in your region.
- If you plan to not use the device for a long time, remove the batteries from the device to avoid the risk of damaging the thermometer.
- Do not insert new and partially used batteries into your device at the same time. This may damage the device.

#### WARNING!

- Keep the thermometer away from children's reach.
- Do not throw batteries into the fire.
- Thermometer can not replace the proper examination and care performed by a medical doctor.

### 02. MEANINGS OF SYMBOLS

The following symbols can be found in the manual, on the device or on the device packaging.

	Applied BF parts
	Do not dispose of the product in the mixed municipal waste container. Dispose of the product in accordance with the guidelines for the disposal of electronic devices of this type
	Follow the instructions for use
	Air moisture range the device properly works in
	<b>Warning</b> The product is resistant to solid objects with a diameter of 12.5mm or more getting inside its casing and against vertically falling drops of water when the device is raised at an angle of up to 15 degrees
	<b>IP22</b> The transport packaging should be dry
	CE mark: The product complies with the requirements of Directive 93/42/EEC
	The transport packaging should not be exposed to sunlight

### 6. CLEANING AND STORAGE

- Do not touch or press the device sensor with your fingers or other objects.
- Use a 75% alcohol swab to wash the thermometer's casing. Do not allow any liquid to get inside the thermometer. Never use corrosive agents, solvents or gasoline to clean the appliance. Do not submerge the product under water or other liquids.
- Keep the thermometer in a dry place, away from dust, dirt, moisture, high temperatures and direct contact with sunlight.

### 7. TYPICAL HUMAN BODY TEMPERATURE

The human body is a complex biological system and temperature range that can be considered "normal" largely depends on what part of the body we measure and on factors such as age, gender, skin color or skin thickness. The thermometer is used to measure body temperature from the forehead of the test subject. The forehead temperature of a healthy person will range from 34.7 °C to 37.8 °C. Women's body temperature is higher than men's by about 0.3°C. In addition, the body temperature of women increases by an additional 0.3-0.5°C during ovulation.

### 8. SPECIFICATION

<b>Measure Mode:</b>	BODY (forehead), SURFACE (surface)
<b>Units of Measure:</b>	Degrees Celsius (°C) or Fahrenheit (°F)
<b>Working temperature:</b>	16-35 °C
<b>Working air humidity:</b>	15-80 %
<b>Storage temperature:</b>	-20 to 55 °C
<b>Air humidity storage:</b>	15-93 %
<b>Measure Distance:</b>	15 to 50 mm from the measurement point
<b>Measure range:</b>	34-43 °C (in BODY mode) 0-100 °C (in SURFACE mode)
<b>Accuracy of measurement:</b>	+/- 0.2 °C (in the range of 35-42 °C); +/- 0.3 °C in other ranges
<b>Memory:</b>	The last 34 measurements taken
<b>Dimensions:</b>	150x37x16.3 mm
<b>Weight:</b>	104 g (without battery)
<b>Batteries:</b>	2x AAA (DC 3V)
<b>Automatic shutdown:</b>	Up to 10 seconds from the last activity

### 9. ELECTROMAGNETIC COMPATIBILITY

Guidance and declarations by the manufacturer on electromagnetic emissions for all equipment and systems.

Guidance and manufacturer's declaration - electromagnetic emissions		
The infrared thermometer shall be adapted for use under the electromagnetic conditions specified below. The thermometer user should make sure that the device is used under these conditions.		
Emission test	Standards met	About the electromagnetic environment
RF emissions CISPR 11	Grupa 1	The infrared thermometer uses RF energy only for its internal functions. Thus, the RF emissions of the device are very low and should not interfere with the operation of nearby electronic devices.
RF emissions CISPR 11	Klasa B	The infrared thermometer can be used in all buildings, including home environments and buildings connected to a public low-voltage power supply network that provides power to buildings for residential use.
Harmonic emissions IEC 61000-3-2	N/A	
Voltage fluctuations/flickering emissions IEC 61000-3-3	N/A	

	Indicates the side of the transport package that should be facing up
	LOT number
	Date of manufacture
	Authorised representative in the European Community
	Manufacturer
	Indicates the temperature range at which the transport packaging should be stored
	Serial number
	Power button

### 03. PRODUCT DESCRIPTION

#### PURPOSE OF THE DEVICE

The infrared thermometer is used to measure the body temperature on the patient's forehead or measurement object temperature. The device is suitable for both home use and use as a medical device. The thermometer can be used to measure temperature regardless of the age of the test subject.

#### DEVICE DESCRIPTION

- Sensor - point towards the person or object being tested.
- Trigger - pressing this button starts the device and performs the measurement.
- "UP" key - used to view saved measurements.
- SET button - Hold the button to turn device sound on/off.
- "DOWN" button - used to view saved measurements.
- MODE button - pressing this button switches the device between the body temperature measurement mode and the surface temperature measurement mode.
- Display - Shows the measurement taken.
- Battery cover - can be removed to replace the batteries in the device.
- Speaker

#### DISPLAY

- Body temperature measurement mode enabled
- Surface temperature measurement mode enabled
- Measured temperature
- Sound on/off
- Low Battery Warning


Guidance and manufacturer's declarations on fault immunity for all equipment and systems.

Guidance and manufacturer's declaration - electromagnetic immunity			
The infrared thermometer shall be adapted for use under the electromagnetic conditions specified below. The thermometer user should make sure that the device is used under these conditions.			
Immunity test	EN 60601 Measured level	Level of compliance	About the electromagnetic environment
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6 kV contact +/- 8 kV air	+/- 6 kV contact +/- 8 kV air	Floors should be made of wood, concrete or ceramic tiles. If the floor is covered with synthetic material, the relative air humidity shall be at least 30%.
IEC 61000-4-4	Power cords: +/- 2kV Input/output cables: +/- 1kV	N/A	The quality of mains power should correspond to a typical hospital or commercial environment.
Surges IEC 61000-4-5	Differential mode: +/- 1kV Common mode: +/- 2kV	N/A	The quality of the mains power should correspond to a typical hospital or commercial environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5% Ut (>95% dip in Ut) For 0.5 cycle  40% Ut (60% dip in Ut) For 5 cycles  70% Ut (30% dip in Ut) for 25 cycles  < 5% Ut (>95% dip in Ut) For 5 seconds	N/A	The quality of the mains power should correspond to a typical hospital or commercial environment. If the user of the infrared thermometer needs to continue using the device when the mains power supply is interrupted, it is recommended that the infrared thermometer has a separate, uninterruptible power source or batteries.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	The frequency of magnetic fields should be at a level corresponding to a typical hospital or commercial environment.
NOTE Ut is the AC mains voltage before applying the test level.			

Recommended separation distances between portable RF communication devices and DEVICE or SYSTEM - for devices and systems that are not life-supporting

Recommended separation distance between portable RF communication devices and infrared thermometer	
The infrared thermometer should be used in electromagnetic environments where RF emission disturbances are controlled. The user of the infrared thermometer can help prevent electromagnetic interference by maintaining a minimum distance between the RF communication devices and the infrared thermometer, in accordance with the following recommendations based on the maximum output power of the communication devices.	Recommended distance between devices depending on transmitter frequency (m)

- Measurement result - normal temperature/ low fever/ high fever
- Degrees Celsius
- Degrees Fahrenheit
- Number of measurements in memory
- Measurement being saved in memory
- Control icon

#### DISPLAY BACKLIGHT

Depending on the measured temperature, the display screen has three colours of backlight:

- Green at ≤37.5°C
- Yellow at 37.6-38.0°C
- Red at >38.1°C

**WARNING!** The colour settings do not indicate anything about the patient's health.

### 04. USING THE DEVICE

#### 1. Install batteries

- Remove the lid of the battery compartment by moving it in the direction indicated by the arrows.
- Insert two AAA batteries inside. Make sure that the battery poles are turned in the right direction.
- Place the cover back in place and press it in the opposite direction than indicated by the arrows to close the battery compartment.

#### WARNING:

- When placing the batteries, make sure that their poles are turned in the correct direction. Incorrectly installed batteries may damage the product.
- If you will not be using the product for an extended period of time, remove the batteries from the product to avoid damaging the device.
- If the batteries or device are showing traces of leakage or mold, stop using them immediately.
- Do not keep batteries near fire or throw them into fire. This can lead to an explosion.
- Do not store batteries in rooms with high temperatures or high moisture.
- To avoid short-circuits, do not keep batteries and metal objects (such as coins or keys) near electrical appliances.

Maximum output power of the transmitter (W)	150 kHz to 80MHz d=(3.5/V1)(Sqrt P)	80MHz to 800MHz d=(3.5/E1)(Sqrt P)	800MHz to 2.5GHz d=(7/E1)(Sqrt P)
0.01	/	0.12	0.23
0.1	/	0.38	0.73
1	/	1.2	2.3
10	/	3.8	7.3
100	/	12	23

In the case of transmitters with the maximum power not mentioned above, the recommended distance d between the transmitter and the infrared thermometer in meters (m) may be calculated on the basis of the formulae listed above. In place of P, insert the maximum power of the transmitter specified by the manufacturer in watts (W).  
**NOTE:** At 80MHz and 800MHz frequencies, adhere to the distances specified in the higher order frequency ranges (e.g. when the device operates at 80MHz, follow the distance specified for frequency range between 80MHz and 800MHz)  
**NOTE:** These guidelines may not be appropriate for all situations. Electromagnetic waves can be influenced by nearby objects, buildings and people.

Guidance and manufacturer's declarations on electromagnetic immunity for devices and systems that are not life-supporting.

Guidance and manufacturer's declaration - electromagnetic immunity			
The infrared thermometer shall be adapted for use under the electromagnetic conditions specified below. The thermometer user should make sure that the device is used under these conditions.			
Immunity test	EN 60601 Measured level	Level of compliance	About the electromagnetic environment
Przewodzone częstotliwości radiowe IEC 61000-4-6	3 Vrms 150 kHz do 80MHz	N/A	Portable RF communication devices should not be closer to any part of the infrared thermometer (including power cords) than the distance specified by the following formulas depending on the frequency of the transmitter:  d=(3.5/V1)(Sqrt P) for 150 kHz to 80MHz d=(3.5/E1)(Sqrt P) for 80MHz to 800MHz d=(7/E1)(Sqrt P) for 800MHz to 2.5GHz  where P is the maximum output power of the transmitter in watts (W) specified by the manufacturer and d is the recommended distance in meters (m). The magnetic field strength for installed RF transmitters determined by the electromagnetic test of their installation (2) site should be less than their level of conformity in each frequency range. (1) Interference may occur near devices marked with this symbol:
Promieniowane fale radiowe IEC 61000-4-3	3 V/m 80 Mhz do 2.5GHz	3 V/m	
<b>NOTE:</b> At 80MHz and 800MHz frequencies, adhere to the distances specified in the higher order frequency ranges (e.g. when the device operates at 80MHz, follow the distance specified for frequency range between 80MHz and 800MHz) <b>NOTE:</b> These guidelines may not be appropriate for all situations. Electromagnetic waves can be influenced by nearby objects, buildings and people.			

### 2. Preparation for measurement

- To ensure the highest possible temperature measurement accuracy, follow these instructions:
  - When you measure the temperature of a person's forehead, select the "BODY" mode on the thermometer.
  - Before measuring, make sure that there is no hair, sweat, cosmetics or clothing on the subject's forehead that could distort the measurement.
  - If the body temperature of the test subject differs significantly from that in the room where the measurement is taken, the test subject should wait at least 5 minutes in the measuring room before starting the measurement.
  - A cold compress or other methods of cooling the forehead for people with fever will make the measured temperature lower than the actual temperature.
  - The temperature in the room where the measurement is performed should be stable. Do not measure in rooms with high airflow, such as rooms cooled with fans or air conditioning.
  - The thermometer should be kept in the room where the measurement is made. If the thermometer has been brought in from another room, leave it in the measuring room for at least 20 minutes before measurement.
  - Do not expose the infrared thermometer to strong sunlight.
  - When measuring, aim the thermometer at the center of the test subject's forehead, above the eyebrows. Keep the thermometer at a distance of ≤ 1cm from the subject's forehead. When you press the "trigger" of the thermometer, the temperature measurement will be shown on the display screen.

### 3. Performing measurement

- Aim the thermometer sensor at the forehead of the person whose temperature you want to measure or at the object whose temperature you want to measure. Press the trigger. The device performs system diagnostics at start-up, so for about two seconds before measuring, the display will look like this:

- After the device diagnostics are done, you will hear a "beep" signal (if the device sound is turned on), indicating that the measurement is correctly performed. The temperature measurement result will be shown on the display.

#### WARNING:

**BODY mode:** If the measured temperature is above 37.6°C, the LCD display turns yellow. If the measured temperature is above 38.1°C, the LCD display turns red and three alarm signals sound from the device. After five measurements taken in quick succession, wait at least 20 seconds before the next measurement to ensure that the measurement is accurate.

#### 4. Change the measurement mode

When the device is on, press the MODE button to change the measurement mode between BODY temperature measurement and SURFACE temperature measurement. The body temperature measurement mode is used to measure the temperature of a human body.

- The strength of magnetic fields of installed transmitters such as radio towers, telephone towers, amateur radio stations, FM/AM waves or television broadcasts cannot be provided at the theoretical level. To assess the electromagnetic environment and the impact of nearby installed RF transmitters, an electromagnetic test must be carried out. If the strength of the electromagnetic fields at the site exceeds the infrared thermometer compliance levels specified above, the thermometer shall be monitored to ensure that it functions normally. If the device is not functioning normally, it may be necessary to move the device.
- In the frequency range from 150kHz to 80MHz, the magnetic field strength should be below 3V/m.

### 10. WARRANTY CARD

The product is covered by a 24-month warranty. The terms of the guarantee can be found at: <https://neno.pl/gwarancja>

Complaints should be reported using the complaint form located at: <https://neno.pl/kontakt>

Specifications and contents of the kit are subject to change without notice. We apologise for any inconvenience.

KGK Trade declares that this device NENO Medic T05 is compliant with significant requirements of directive 2014/53/EU. The text of this declaration may be found on the website: <https://neno.pl/download/DOC/deklaracja-CE-Neno-Medic-T05.pdf>

The crossed out trash can symbol indicates that unusable electrical or electronic devices, its accessories (such as power supplies, cords) or components (for example batteries, if included) cannot be disposed of alongside with household waste. In order to dispose of the devices or its components (for example, batteries) deliver the device to the collection point, where it will be accepted free of charge. Disposal is subject to the latest version of the WEEE Directive (2012/19 / EU) and the Directive on batteries and accumulators (2006/66 / EC). Proper disposal of the device prevents degradation of the natural environment. Information about the collection points of the facilities is issued by the competent local authorities. Incorrect disposal of waste is subject to penalties provided for by the law in force in the given area.

**Importer:**  
KGK Trade sp. z o.o. sp. k.  
Os. Urocz 12  
31-953 Cracow  
Poland

**EU Authorized Representative:**  
 Prolix GmbH,  
Brehmstr. 56, 40239, Düsseldorf,  
Germany

**Manufacturer:**  
 Xuzhou Yongkang Electronic Science Technology Co., Ltd.  
4F building C8, 40, Jingshan Road, Economic and  
Technological Development Zone,  
221000 Xuzhou, PEOPLE'S REPUBLIC OF China

BY **Xblitz**.

**neno**  
www.neno.pl  
Designed by Futurus

The quality of the product is guaranteed by the  
TUV SUD Product Service GmbH

**CE 0123**

Last update of the user's manual: 24.09.2020