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## sente

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
**Wytwórcza:**  
Alic Medical Shenzhen, Inc  
4/F, B Building,  
ShenJubao Modern Optical Factory,  
Kangxi Street Pingshan District, 518122 Shenzhen City,  
People's Republic of China

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BY **Xblitz.**

**neno**  
www.neno.pl  
inspired by Children,  
designed by Parents

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# nenosente

## USER'S MANUAL


Dear Customer,  
thank you for purchasing your Neno Sente. The compressor nebulizer you have purchased is a medical device used to administer medicine assigned or prescribed by your doctor to the respiratory tract. The nebulizer is powered by an alternating current of 220V. The device should be used under the supervision of a doctor or therapist.  
Before using the product, please read the instructions below carefully.

### 1. SAFETY MEASURES

- Before use, carefully read the user manual and keep it for future reference.
- The device may only be used for the purposes specified in the manual. The manufacturer is not responsible for damages caused by improper use.
- The product is not suitable for use on infants and newborns.
- The nebulizer is NOT used for the nebulization of oil-based medicine.
- The device housing cannot be opened.
- The device is not suitable for the administration of anesthetics and must not be part of the lung ventilator system.
- The device is not to be used for the administration of drugs in the form of a suspension or drugs with high viscosity. Contact the pharmacy that sold you that medicine for information on proper administration.
- Make sure that the device is not visibly damaged before each use. If the device is damaged, stop using the device and contact the service.
- Before each use of the nebulizer, check the device, air filter, mask and rest of the equipment. Dirty or worn parts should be replaced with new ones.
- The patient's face mask is treated as plastic for the purposes of disposal after washing and disinfecting it.
- Do not disassemble the product. This will void the warranty.
- Keep the device and its packaging out of the reach of children.
- If the device plug does not fit into the power outlet, stop using the device and contact service. Do not use adapters or power splitters.
- Unplug the device if it is not used or if it is left unattended.
- The device complies with applicable EU directives and has the CE marking.
- If the nebulizer is used frequently, it is recommended to replace the device every 6 months. Replace the device earlier if it clogs.
- The device must not be used by persons with limited physical, mental or cognitive abilities, or who do not have the appropriate knowledge or experience to use the nebulizer correctly (including children). In order for such persons to be able to use the device, they must do so under the supervision of the person responsible for them that has read the user's manual.
- Keep the power cord of the device away from hot surfaces.
- Use only accessories recommended by the manufacturer.
- Keep your device away from children. They can swallow the small elements that make up the device set.
- The device is not waterproof. Do not wash it under running water or submerge it in liquids. Do not touch the device with wet hands.
- If the device falls into the water, immediately disconnect it from the power supply.
- Do not use the device in rooms containing flammable gases, compressed oxygen or other flammable spray products.
- The spray of the medicine should be kept away from eyes.
- The capacity of the medicine container is up to 8ml, this dose must not be exceeded.
- Do not block the ventilation openings of the device. Do not place the nebulizer on soft surfaces that may block the ventilation openings.
- Before cleaning, unplug the device from the power supply.
- Always remove the unused amount of medicine from the device after use. For each inhalation, use a fresh dose of medicine.
- Store your device and its accessories in a dry place. Moisture can lead to the development of bacteria.
- Keep the device away from water. Do not use the device while bathing.
- Do not expose the device to direct sunlight or high temperatures.
- Each patient should have a personal set consisting of a nebulizer, mouthpiece and mask. This allows you to maintain the highest level of hygiene for all patients treated with the device.

### 2. SET ELEMENTS

### 3. PACKAGE CONTENT

- Compressor
- Nebulizer
- 2 masks (one for an adult and one for a child)
- Mouthpiece
- Air tube
- Five pcs. interchangeable filters
- User's manual

### 4. USING THE DEVICE

- The accessories included are sterile and are used for a single sterile inhalation. Wash the accessories thoroughly with water before using them next time. Each patient should have a personal nebulizer, mouthpiece and mask.
- Place the compressor nebulizer on a flat, stable surface. Make sure you can easily reach your device while you're sitting.
  - Remove accessories from the bag
  - NOTE: Before first use, carefully refer to the instructions.
  - Gently unscrew the top of the nebulizer counterclockwise to take it apart.
  - Make sure there is a plastic cone at the bottom of the container. Pour the assigned dose of the drug into the container.
  - Place the top of the nebulizer in the container. Make sure the parts fit well together. Gently screw the nebulizer parts together by turning the upper part of the nebulizer clockwise.
  - You can touch the power switch and filter cover if the contact time is less than 10 seconds.
  - Connect one of the ends of the tube to the nebulizer.
  - Connect the other end of the tube to the designated tube connector between the power button and the air filter.
  - Connect the mouthpiece and the face mask of your chosen size to the nebulizer.
  - Plug the device into an electrical outlet. Make sure your device is OFF before you do.
  - Press the power switch to turn on the device and start therapy.
  - After the therapy is complete, turn off the device and disconnect it from the power outlet.

NOTE: The compressor engine is equipped with thermal protection that turns off the device before it can overheat. If the device has switched off on its own:  
a. Turn off the device.  
b. Unplug the device.  
c. Wait 30 minutes before next use. Make sure that the ventilation openings of the unit are not blocked.

If the device does not start, make sure it is plugged in. If the device continues not working after connecting it to the power supply, contact the authorized service.

### 5. CLEANING

- Compressor cleaning:
  - Wipe the device with damp cloth every day.
  - Do not use powdered cleaners or soaps as they may damage the finish of the machine.
- Cleaning the tube, nebulizer, mouthpiece and face mask:
  - Wash the components with running water after using the appliance.
  - Remove the remnants of the medicine from the accessories.
  - Let the accessories dry naturally. Make sure they are dry before storage – moisture can lead to the growth of bacteria in the device.
- Filter replacement:
  - Do not use cotton or similar materials. Do not wash the filter. Do not use the device without a filter.
  - Replace the air filter every 60 days or when the filter gets dirty or wet.
  - To replace the filter:
    - remove the filter cover,
    - replace the used filter with a new one,
    - insert the filter cover back in place.

### 6. USED ICONS

The following markings may be included in the instructions, on the packaging, on the device or on the accessories for the device. Some symbols represent standards and compliance for the device and its use.

BF-type parts used

SN Device serial number

Manufacturer

Official representative in the European Community

CE mark: The product complies with 43/42/EEC directive

Disposal: Do not dispose of this product in a mixed municipal waste container. The product must be disposed of in accordance with the regulations in your region.

Follow the instructions for use

IP21 The product is resistant to water and particle matter getting inside its casing

### 7. ELECTROMAGNETIC COMPATIBILITY

The ME DEVICE is adapted for use in home and similar environments for healthcare purposes.  
**Warning:** Do not start the device near high-frequency surgical devices or in magnetic resonance rooms with RF shielding, as there are high amounts of electromagnetic disturbances in these locations.

**Note:** Using the device near other equipment may lead to measurement irregularities and should be avoided. If the use of two devices in each other's vicinity is unavoidable, both devices should be carefully monitored to ensure that they are working properly.

If any: a list of all cables and maximum cable length (if applicable), transducers and other ACCESSORIES which are replaceable by THE RESPONSIBLE ORGANIZATION and which may affect the compliance with the standards of ME EQUIPMENT or ME SYSTEM with the requirements set out in Clause 7 (EMISSIONS) and Clause 8 (IMMUNITY). ACCESSORIES can be specified in general (e.g. shielded cable) or specifically (e.g. specifying the MANUFACTURER and TYPE OF EQUIPMENT).

If any: the functioning of medical equipment or medical system has been defined as ESSENTIAL PERFORMANCE and a description of the changes that the OPERATOR may experience if ESSENTIAL PERFORMANCE is disturbed by EM DISTURBANCES (the term ESSENTIAL PERFORMANCE does not need to be used).

- All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to EM disturbances for the period of the device's expected service life.
- Guidance and manufacturer's declaration – electromagnetic emissions and immunity.

Table 1

Guidance and manufacturer's declaration – electromagnetic emission	
Emission test	Standard met
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	N/A
Voltage fluctuations/flicker emissions IEC 61000-3-3	N/A

Table 2

Guidance and manufacturer's declaration – electromagnetic immunity		
Immunity test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact +/- 2 kV +/- 4 kV, +/- 8kV, +/- 15kV air	+/- 8 kV contact +/- 2 kV +/- 4 kV, +/- 8kV, +/- 15kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply cables: +/- 2kV Input/output cables: +/- 1kV	N/A
Surges IEC 61000-4-5	Line(s) to line(s): +/- 1kV Line(s) to earth: +/- 2kV Repetition frequency 100 kHz	N/A
Voltage dips, short interruptions and voltage variations in power supply cables IEC 61000-4-11	0% 0.5cycle At 0°, 45°, 90°, 135°, 180°, 225°, 275° and 315° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycles	N/A
Magnetic field frequency IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted radio frequencies IEC 61000-4-6	150kHz to 80MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	N/A
Radiated radio waves IEC 61000-4-3	10 V/m80MHz – 2,7GHz±0% AM at 1kHz	10 V/m 80MHz – 2,7GHz 80% AM at 1kHz

NOTE: Ut is the AC mains voltage before applying the test level.

Table 3

Guidance and manufacturer's declaration – electromagnetic immunity						
Test frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (on/off)	Distance (m)	Measured Immunity level (V/m)
385	380-390	TETRA 400	Pulse modulation 18Hz	1.8	0.3	27
450	380-390	GMRS 460, FRS 460	FM +/- 5 kHz, deviation 1 kHz sine	2	0.3	28
710 745 780	704-787	LTE band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9
810 870 930	800-960	GSM 800/900, TETRA 800, IDEN 820, CDMA 250, LTE Band 5	Pulse modulation 18Hz	2	0.3	28
1720 1845 1970	1700-1990	GSM 1900; CDMA 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217Hz	2	0.3	28
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Psmo 7	Pulse modulation 217Hz	2	0.3	28
5240 5240 5785	5100-5800	WLAN 802.11 a/h	Pulse modulation 217Hz	0.2	0.3	9

Radiated RF IEC 61000-4-3 (test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communication equipment).

### 8. TECHNICAL SPECIFICATION

Power supply: AC 220 V 50 Hz  
Power consumption: 180W A  
Capacity of the medicine container: 2 ~ 8 ml  
Particle size: 1 up to 5 um  
MMAD: 3.2 um  
Noise level: <= 70 dBA  
Amount of nebulized fluid : >= 0.15 ml / min  
Operating pressure range: 80 up to 120 kPa  
Air flow: >= 5 l / min  
Operating temperature: 10 - 40 °C  
Operating atmospheric pressure: 860 - 1060 hPa  
Storage temperature: -20 up to 55 °C  
Storage atmospheric pressure: 500 - 1060 hPa  
Storage/operating air humidity: 15 ~ 93 % RH  
Dimensions (compressor): 170x150x87.5 mm  
Weight (compressor): 1150 g  
Weight (full set): 1240 g  
Air tube length: 1538 m

### 9. WARRANTY CARD

The product is covered by a 24-month warranty. The terms of the guarantee can be found at: <https://neno.pl/gwarancja>

Complaints should be reported using the complaint form located at: <https://neno.pl/kontakt>

Specifications and contents of the kit are subject to change without notice. We apologise for any inconvenience.

KGK Trade declares that this device Neno Sente is compliant with significant requirements of directive 2014/53/EU. The text of this declaration may be found on the website: <https://neno.pl/download/DOC/deklaracja-CE-Neno-Sente.pdf>

The crossed out trash can symbol indicates that unusable electrical or electronic devices, its accessories (such as power supplies, cords) or components (for example batteries, if included) cannot be disposed of alongside with household waste. In order to dispose of the devices or its components (for example, batteries) deliver the device to the collection point, where it will be accepted free of charge. Disposal is subject to the recast version of the WEEE Directive (2012/19 / EU) and the Directive on batteries and accumulators (2006/66 / EC). Proper disposal of the device prevents degradation of the natural environment. Information about the collection points of the facilities is issued by the competent local authorities. Incorrect disposal of waste is subject to penalties provided for by the law in force in the given area.

EU Authorized Representative:  
Shanghai International Holding Corp. GmbH (Germany)  
Address: Eiffelstrasse 80, 20537 Hamburg, Germany

Manufacturer:  
Alican Medical Shenzhen, Inc  
4/F, B Building,  
Shentubao Modern Optical Factory,  
Kengzi Street Pingshan District, 518122 Shenzhen City,  
People's Republic of China

The quality of the product is guaranteed by the  
TUV SUD Product Service GmbH

Last update of the user's manual: 21.08.2020

Importer:  
KGK Trade sp. z o.o. sp. k.  
Os. Urocz. 12,  
31-953 Cracow  
Poland

Made in P.R.C.

BY Xblitz.  
  
www.neno.pl  
Designed by Children,  
Designed by Parents

# nenosente